



Chair Glide for ADA Compliance









Chair Glide with Tapered Track



Heavy-duty stops on base plate limit distance of travel.

Tapered track improves access to patient.

Locking pedal secures chair in any position along its 27" path.

Case-hardened steel rollers allow you to easily slide chair out of the way for wheelchair access.



The high-density polyethylene plate fits on any floor or carpet and assures quite travel. Centered on the plate, a track guides the carriage and prevents lateral movement.

Carriage Fits All Reliance Chairs

A heavy-duty, cast aluminum carriage fits any Reliance exam chair.

Instant ADA Compliance

Along with a low-base Reliance chair, the 940 chair glide will help you comply with the mandatory requirements of the Americans with Disabilities Act.

Potential Tax Advantages

Apply the entire cost of the 940 to the tax credit provided by the ADA legislation. Check with your tax professional for more details.

Tech Specs

Base plate length: 52" overall; wide section: 28"; narrow: 24". Base plate width: wide section: 20"; narrow: 6". Track: ¾" wide by 3/16" high by 52" long. Overall carriage length: 29"; width: 26".



